

Dual Oxygen Cart

HCPCS:




FIGURE A



- Lightweight silver vein cart can accommodate two "D" or "E" size cylinders
- Handle is adjustable to accommodate most individuals
- Silver vein finish easy-to-maintain & less institutional looking (Figure A)
- Large, composite wheels allow for easy, safe maneuvering

PRODUCTS

Description

13001SV-2	Oxygen CylinderCart,Dual,Silver,2/cs
-----------	--------------------------------------

SPECIFICATIONS

Description

Warranty	Limited Lifetime
Weight Capacity	80 lbs
Height	12.2"
Length	38.6"-41.9"
Carton Shipping Weight	22 lbs.
Weight	5.4 lbs.