

Meridian Ultra-Care 4800 8" Mattress with High-Output Pump

HCPCS:

- The Ultra-Care 4800 is a combination mattress system providing both alternating pressure and low air loss for the management and treatment of pressure ulcers.
- Quick disconnect cells allow any area of the mattress to deflate and be zero pressure without affecting alternating benefits.
- 8" alternating pressure mattress with low air loss regardless of operating mode
- Pillow function maintains the air cells at the head for patient comfort
- Air cells are made with nylon/polyurethane for increased durability and reliability and can quickly disconnect for easy replacement
- CPR quick deflate at head of mattress
- Quilted nylon zipper cover is low shear, fluid-resistant, and vapor permeable
- High-output digital pump has an easy-to-use design and offers both static and alternating modes of therapy for maximum comfort
- Includes 8" mattress, cover, and high-output pump
- Weight capacity: 450 lbs
- Dimensions: 80" x 36" x 8"

PRODUCTS

Description

| | |
|-------|----------------------------------|
| 4800E | ULTRA-CARE 4800E APM SYSTEM 1EA/ |
|-------|----------------------------------|
