

Low Bed Overbed Table

HCPCS: E0274



- Designed specifically for use with low beds with 3.25" frame to floor clearance
- Walnut woodgrain melamine top
- Top can be raised or lowered in infinite settings between 19.5"- 27.5"
- Tabletop can be raised with the slightest upward pressure
- The frame is constructed with tubular steel and durable, attractive chrome finish
- Swivel casters allow for easy maneuverability

| PRODUCTS | Description |
|----------|------------------------------|
| 13081 | Overbed Table,Low Style,1/ea |

| SPECIFICATIONS | Description |
|------------------------|------------------|
| Warranty | Limited Lifetime |
| Weight Capacity | 100 lbs |
| Height | 19.5" - 27.5" |
| Base Dimensions | 15" x 28" |
| Carton Shipping Weight | 28 lbs. |
| Weight | 25 lbs. |