

## Pivot & Tilt Overbed Table

HCPCS: E0274



FIGURE A

- The mast of the table pivots and can be locked in one of the three positions from flat on the floor to 90° (Figure A)
- The tabletop tilts with the angle of the mast to provide a flat surface regardless of the angle of the mast
- The pivot feature allows the table top to be positioned closer to individuals in a bed, sitting in a wheelchair, or folded for storage
- Four 2" casters (two locking)
- Only wheels have to be installed
- New and improved larger tubing

PRODUCTS	Description
13000	Overbed Table,Pivot/Tilt,1/cs
SPECIFICATIONS	Description
Warranty	Limited Lifetime
Weight Capacity	40 lb
Carton Shipping Weight	22 lb
Tabletop Adjustment	28" - 37"
Weight	21 lb