

PB42BARBED Foot spring, electronics, pendant, crank, side support bars

HCPCS:

- The Bariatric Bed Foot Springs allow bariatric patients to adjust their foot position to their liking.
- For use with ProBasics 42" Bariatric Bed (#PB42BARBED).
- Foot spring section for ProBasics Bariatric Bed (#PB42BARBED).
- To be used with Bariatric Head Spring section (#PB3000) and Bariatric Bed Ends (#PB3301).
- Motorized positioning of foot section provides maximum comfort and versatility to patient and caregiver
- Easy-to-use hand control allows for multiple bed positioning for patients
- Easy set up; no tools necessary
- Weight capacity: 600 lb
- For a wide range of foot position options for bariatric patients, purchase the Bariatric Bed Foot Spring.

PRODUCTS

Description

PB3010	BARIATRIC BED, 42", FOOT SPRNG 1EA/
--------	-------------------------------------