

ProBasics Bariatric Trapeze (Complete Unit)

HCPCS:



- The ProBasics Bariatric Trapeze features easy assembly and a lightweight design.
- Quick onsite set-up and the maximum user weight capacity of 650 lbs., makes this the best choice for larger users.
- This trapeze is designed to provide support, increase stability and assist the user when re-positioning himself/herself in bed.
- To allow for maximum versatility, the overhead triangle grip with adjustable chain can be positioned in several locations along the overhead boom.
- FEATURES
- Helps users to change positions while in bed.
- Overhead triangle grip and adjustable chain can be positioned in several locations along the overhead boom - just slide the slide ring to position.
- Easy-to-assemble with the required tools.
- Drill, rubber mallet and socket with wrench are required.
- The overhead boom swivels from 0° to 180° (for storage only)
- Maximum user weight capacity of 650 lbs.

PRODUCTS	Description
PBBT	BAR TRAPEZE - 650 LB WT CAP 1EA/