

ProBasics Bed Cart for 36" wide Homecare Beds

HCPCS:

- The ProBasics Bed Cart has been designed for the efficient coordination of homecare bed fleets.
- This cart ships fully assembled and uses the homecare bed ends to create an easy-to-move and compact storage system.
- The strong, sturdy shelf holds the springs, hardware and user guide for a quick turn-around.
- Only sold in case packs of 2 units.
- Steel frame fits onto most standard (36 wide) bed ends
- Composite shelf is rugged for years of dependable service
- Keeps your warehouse organized
- Technicians can roll the beds from warehouse-to-truck-to-customer
- Helps to manage cleaning routines and sanitation protocols

PRODUCTS

Description

PBHBCART36

HEMOCARE BED CART W/ PLATFORM 2 EA/CV