

ProBasics Bed Mount Trapeze

HCPCS:



- The ProBasics Bed Mount Trapeze provides weaker patients with something secure to grasp and lift themselves up with when changing positions in bed.
- It has a strong, sturdy base that can be mounted directly to the home care bed with the included bed clamps.
- Attractive, chip resistant finish on tubular steel
- Adjustable height, horizontal position and hand bar support distance to assist a range of patients
- Attaches easily and securely to headboard
- Easy, tool-free installation
- Weight capacity: 250 lbs
- Includes trapeze and clamps for bed attachment
- Also available in free standing design (order item PBST)
- To help patients independently move in and out of bed, purchase the ProBasics Bed Trapeze.

PRODUCTS	Description
PBST-TRPZ	BED-MOUNT STD TRAPEZE W/CLAMPS 1EA/