

BESTCARE PL273 CONTROL BOX 1EA/

HCPCS:

PRODUCTS	Description
DSC700	BESTCARE DIGITAL PATIENT 1EA/
KIT-TA23-PL400	BESTCARE UPGRADE KIT-PL400 TO ELECTRIC
PBPLSM001	PATIENT LIFT HYDRAULIC 1EA/
PBST	FLOOR MOUNT TRAPEZE 1EA/
PL182	BESTCARE DELUXE ELECTRIC PATIENT LIFT
PL228	BESTCARE PREMIUM ELECTRIC PATIENT LIFT
PL400EF	BESTCARE ELECTRIC FOLDING PATIENT LIFT
PL400H	BESTCARE HYDRAULIC PATIENT LIFT
PL400HE	BESTCARE GENESIS ELECTRIC PATIENT LIFT
PL500	BESTCARE PATIENT LIFT
SA182	BESTCARE DLX ELECTRIC STAND ASSIST LIFT