

ProBasics Four Bar Clamp-On Half Rail

HCPCS:



- The ProBasics Four Bar Clamp-On Half Rail is perfect for patients who will be getting in and out of bed by themselves, as they will be able to reach the release to drop the rails when exiting bed.
- It is also recommended for patients that don't move around in bed too much.
- No gap style reduces entrapment hazards and adds protection
- Easy to install; no tools required for installment
- Fits all woven and spring style hospital beds
- Durable and constructed of welded steel with brown vein finish
- Dimensions: 19"H x 37"L
- For added patient protection, purchase your ProBasics Four Bar Clamp-On Half Rail.

PRODUCTS	Description
PB7035	HALF LENGTH BED RAILS (PR) 1PR/