

ProBasics IV Pole, 4-Leg, 2-Hook

HCPCS:

- The ProBasics IV Pole is perfect for hospital or at-home IV therapy, with twin wheel swivel casters that are easy to glide across hard floor or carpet.
- It is also accessible for everyone, with a wide height adjustment range.
- Chrome-plated and made of durable hammertone steel
- Offset low profile base keeps IV pole stable
- Twin wheel swivel casters allow for IV pole to move smoothly across hard floor or carpet
- Twist lock provides full range of height positions
- Two-pronged removable top keeps IV containers securely in place while making IV pole easy to store and transport
- Height adjustment: 45.5" to 79"
- Base size: 23"
- Quantity: 2/case
- Purchase your ProBasics IV Pole for smooth mobility during hospital or at-home IV therapy.

PRODUCTS	Description
PB5035	IV POLE, 4 LEG, 2 HOOK 2EA/CV