

## Automatic Blood Pressure Monitor, Wrist Model

HCPCS: A4670-1,A4670-1



FIGURE A

- Large easy-to-read digital display
- 120 memory recall
- Date & time display
- Pulse Rate Indicator
- Auto-off function
- Lightweight & easy to use
- Quiet inflation method
- Wrist model

### PRODUCTS

### Description

BP2116	BP Monitor Wrist Auto 1 EA
--------	----------------------------

### SPECIFICATIONS

### Description

Warranty	2 Years
Measurement Range	Pressure: 0~300mmHg +/-3mmHg
Power Supply	2 x AAA Batteries
CA Prop 65	[WARNING: This product can expose you to chemicals that are known to the State of California to cause cancer and birth defects or other reproductive harm. For more information go to <a href="http://www.P65Warnings.ca.gov">www.P65Warnings.ca.gov</a> .]