

Deluxe Automatic Blood Pressure Monitor, Upper Arm

HCPCS: A4670-1,A4670-1



FIGURE A



FIGURE B

- WHO Classification
- Irregular Heart Beat Detection
- One Touch Operation
- X-Large Easy-to-Read Screen
- Date & Time Display
- Pulse Rate Indicator
- Digital Error Messages
- Travel Bag
- Auto-Off Function
- Memory stores 120 readings (2 x 60)
- Average Last 3 Results

PRODUCTS

Description

BP3400	BP Monitor, Auto, Upper Arm, 1/EA
--------	-----------------------------------

SPECIFICATIONS

Description

Warranty	2 Years
Measurement Range	Pressure: 0mmHg~300 mmHg Pressure: +/- 3 mmHg
Power Supply	4x AA Batteries or AC Adapter
CA Prop 65	[WARNING: This product can expose you to chemicals that are known to the State of California to cause cancer and birth defects or other reproductive harm. For more information go to www.P65Warnings.ca.gov .]

ACCESSORIES

Description

BP2001	Cuff,Std,BP3400,BP3600,BPHL,8.5-14, EA
BP2610	A/C Adapter, BP3400,BP3600,1/ea, EA