

## Deluxe Automatic Blood Pressure Monitor, Wrist

HCPCS: A4670-1,A4670-1



FIGURE A



FIGURE B

- WHO Classification
- Irregular Heart Beat Detection
- One Touch Operation
- Large Easy-to-Read Screen
- Date & Time Display
- Pulse Rate Indicator
- Digital Error Messages
- Travel Bag
- Auto-Off Function
- Memory Stores 120 Readings (2 x 60)
- Average Last 3 Results

### PRODUCTS

### Description

BP3200	BP Monitor, Deluxe, Auto,Wrist, 1/EA
--------	--------------------------------------

### SPECIFICATIONS

### Description

Warranty	2 Years
Measurement Range	Pressure: 0mmHg~300mmHg Pressure: +/- 3mmHg
Power Supply	2x AAA Batteries

### ACCESSORIES

### Description

BP2610	A/C Adapter, BP3400,BP3600,1/ea, EA
BP2000	Cuff,PlusSize,12.5"x17.75",BP3400,3600, EA
BP2001	Cuff,Std,BP3400,BP3600,BPHL,8.5-14, EA