

Deluxe Automatic Blood Pressure Monitor, Wrist

HCPCS: A4670-1,A4670-1



FIGURE A



FIGURE B

- WHO Classification
- Irregular Heart Beat Detection
- One Touch Operation
- Large Easy-to-Read Screen
- Date & Time Display
- Pulse Rate Indicator
- Digital Error Messages
- Travel Bag
- Auto-Off Function
- Memory Stores 120 Readings (2 x 60)
- Average Last 3 Results

PRODUCTS	Description
BP3200	BP Monitor, Deluxe, Auto,Wrist, 1/EA
SPECIFICATIONS	Description
Warranty	2 Years
Measurement Range	Pressure: 0mmHg~300mmHg Pressure: +/- 3mmHg
Power Supply	2x AAA Batteries
CA Prop 65	[WARNING: This product can expose you to chemicals that are known to the State of California to cause cancer and birth defects or other reproductive harm. For more information go to www.P65Warnings.ca.gov .]
ACCESSORIES	Description
BP2610	A/C Adapter, BP3400,BP3600,1/ea, EA
BP2000	Cuff,PlusSize,12.5"x17.75",BP3400,3600, EA
BP2001	Cuff,Std,BP3400,BP3600,BPHL,8.5-14, EA