

Economy Automatic Blood Pressure Monitor, Upper Arm

HCPCS: A4670-1,A4670-1



FIGURE A



FIGURE B

- WHO Classification
- Irregular Heart Beat Detection
- One Touch Operation
- Large Easy-to-Read Screen
- Date & Time Display
- Pulse Rate Indicator
- Digital Error Messages
- Travel Bag
- Auto-Off Function
- Memory Stores 120 Readings
- Average Last 3 Results

PRODUCTS

Description

| | |
|--------|--------------------------------|
| BP3600 | BP Monitor,Auto,Upper Arm,1/EA |
|--------|--------------------------------|

SPECIFICATIONS

Description

| | |
|-------------------|---|
| Warranty | 2 Years |
| Measurement Range | Pressure: 0mmHg~300 mmHg Pressure: +/- 3 mmHg |
| Power Supply | 4x AAA Batteries or AC Adapter |
| CA Prop 65 | [WARNING: This product can expose you to chemicals that are known to the State of California to cause cancer and birth defects or other reproductive harm. For more information go to www.P65Warnings.ca.gov .] |

ACCESSORIES

Description

| | |
|--------|--|
| BP2610 | A/C Adapter, BP3400,BP3600,1/ea, EA |
| BP2000 | Cuff,PlusSize,12.5"x17.75",BP3400,3600, EA |
| BP2001 | Cuff,Std,BP3400,BP3600,BPHL,8.5-14, EA |