

Pure Expressions Dual Channel Electric Breast Pumps

HCPCS: E0603-1,E0603-1



FIGURE A



FIGURE B



FIGURE C



FIGURE D



FIGURE E

- LCD Display provides status of programs and functions
- Independent & custom suction for each breast 9 adjustable suction strengths
- 2 Phase Expression Stimulation and Expression Mode
- Includes 2 Soft Silicone Cushion Inserts for Optimum Comfort
- Easy to assemble and clean Silent for discreet pumping
- Closed Loop Anti-backflow System
- Single or Double Pumping Options
- BPA and Latex Free
- Two 4oz Bottle Feeding Adapter Kits and Storage Containers with Lids
- Convenient Carry Bag with Zipper
- Insulated Cooler Bag comes with RTLBP2200 only
- Includes A/C Power Supply & External battery pack 4 AA alkaline batteries (not included)
- One Year Limited Warranty

PRODUCTS

Description

| | |
|-----------|---------------------------------------|
| RTLBP2000 | Double Electric Plus Breast Pump 1/ea |
|-----------|---------------------------------------|

SPECIFICATIONS

Description

| | |
|-----------------------|--|
| Warranty | One Year Limited |
| Product Weight | .9 lbs |
| Product Dimensions | 4" (W) x 6.25" (L) x 2.75" (H) |
| Batteries (pair) | 4 AA Alkaline |
| Contains Latex | No |
| Retail Box Dimensions | RGLBP200 12.2" (W) x 7.5" (L) x 11" (H) RTLBP2200 12.2" (W) x 9.1" (L) x 12.8" (H) |
| Retail Box Weight | 5.4 lbs. |

ACCESSORIES

Description

| | |
|-------|--------------------------------|
| BP001 | Breast Pump Carry Bag,1/ea, EA |
|-------|--------------------------------|