

Pure Expressions Single Channel Electric Breast Pump

HCPCS: E0603-1,E0603-1



FIGURE A



FIGURE B



FIGURE C



FIGURE D

- LCD Display provides status of programs and functions
- 9 Adjustable suction strengths
- 2 Phase Expression Stimulation and Expression Mode
- Easy to assemble and clean silent for discreet pumping
- Closed Loop Anti-backflow System
- Single Pumping
- BPA and Latex Free
- One 4oz Bottle Feeding Adapter Kit and Storage Container with Lid
- Convenient Drawstring Bag
- A/C Power Supply Included Internal battery power source 4 AA alkaline batteries (not included)
- One Year Limited Warranty

PRODUCTS	Description
RTLBP1000	Single Electric Breast Pump 1/ea

SPECIFICATIONS	Description
Warranty	One Year Limited
Product Weight	.5 lbs.
Product Dimensions	4" (W) x 4.75" (L) x 2.25" (H)
Batteries (pair)	4 AA Alkaline
Contains Latex	No
Retail Box Dimensions	8.8" (W) x 6" (L) x 6.9" (H)
Retail Box Weight	2.2 lbs.

ACCESSORIES	Description
BP001	Breast Pump Carry Bag,1/ea, EA
BP0010	Power Adapter, 6V, RTLBP1000, 1/ea, EA