

CAM Walker Air Boot, Small

HCPCS:

- This CAM Walker Air boot offers support, mobility, and protection for any activities during the day, while recovering from a variety of foot and ankle injuries.
- You can use the air pump to contour and adjust to your fit.
- Patients of this boot usually have sprained ankle, stress fracture, postoperative use, metatarsal fracture, Lisfranc injuries, broken foot, Achilles tendonitis.
- FEATURES:
 - Ankle and leg protection
 - Decrease pain symptoms
 - Improve healing time
 - Immobilization of fractures
 - Improve ambulation during post-op and post-surgery
 - Sturdy hook and loop straps in place
 - Skid-resistant 360° rocker sole

PRODUCTS

Description

CAMAIRSM

WALKING BOOT LOW TOP AIR SM 18EA/CV
