

## StrengthTape Kinesiology Tape, 5M, Black

HCPCS:



- STRENGTHTAPE® features long-lasting durability and effectiveness on common injuries.
- Class exclusive AllSport adhesive helps it sticks for up to 7 days through water, sweat, and mud to support injured muscles and tissues while still allowing full range-of-motion.
- STRENGTHTAPE® is made of breathable cotton and Spandex.
- It is not made with natural rubber latex.
- STRENGTHTAPE® Pre-Cut rolls include tape for up to 10 complete applications.
- Uncut rolls come with 5 meters of uncut tape.
- Each roll includes easy-to-follow instructions for the most common injuries.
- Additional application videos can be found at [www.strengthtape.com](http://www.strengthtape.com).

### PRODUCTS

### Description

6305-510PC

ST 5M Precut Roll - Black 144EA/CV