

Heating Pad

HCPCS:



FIGURE A



FIGURE B




FIGURE C



FIGURE D



FIGURE E



FIGURE F



FIGURE G



FIGURE H

- This digital heating pad provides comfortable and fast-acting dry or moist heat within seconds. Featuring six unique heat and time settings, the user can adjust the pad to their individual needs for maximum results. The heating pad is lightweight, flexible and soft to conform to any part of the body.
- Safe, Soothing & Fast-Acting Heat
- Absorbent sponge in pad provides moist heat that penetrates tissue quickly and safely
- Pad heats in seconds and monitors real-time temperature to prevent overheating
- Six temperature settings (113/122/131/140/149/158) and countdown timer settings (0:15/0:30/0:45 /1:00/1:30/2:00)
- Digital display features on/off switch, temperature and timer control, and indicator light
- Machine-washable pad
- Weighs 1.5 lbs. with large dimensions of 12 x 24 x 0.4
- Comes with portable tote bag for storage or convenient on-the-go use and elastic fastening strap

PRODUCTS	Description
RTL19G007GR	Digital Heating Pad, 1EA, 6/CV
SPECIFICATIONS	Description
Warranty	Limited Lifetime
Color	Gray