

CareTec IV

HCPCS:

- The Care Tec IV combines the waveforms of Russian, EMS, TENS and IF therapy in a portable unit.
- The device offers preset modes that can be user defined.
- Maximum portability and simple operation offering you all the options of a clinical device in a portable package.
- Dual channel device
- Simple to operate, portable device for pain management
- Four waveforms, Russian, EMS, TENS and IF
- Comes complete with carrying case, lead wires, 9-volt battery, AC adaptor, manual and pack of four self-adhesive reusable electrodes
- Power Source: Standard 9-volt battery or AC adaptor
- Dimensions: 4.5" x 2.55" x 0.9"
- Weight: 0.28 lb.
- with battery

PRODUCTS

Description

DT4444	CareTec IV - Digital 32EA/CV
--------	------------------------------
