

## EMS 5.0 Muscle Stim Analog Unit

HCPCS:

- NOTE: Class II - Rx Required
- Qty: 32 per master carton Contains: 4/pk E1F1515TC2, 1-9 volt battery TA9050-I, 2/pk WW3005 leads
- The EMS 5.0 Muscle Stim Analog Unit is designed for personal use, however it is strong enough to be used in an office setting.
- The EMS unit uses electrical impulses to contract and relax muscles helping those who experience muscle tension in the upper back and neck areas.
- Can be used to regain muscle strength after surgery or being in a cast.
- The EMS unit is portable and may be clipped to a belt, shirt, pocket or other clothing.
- System includes device, lead wires, four self-adhesive reusable electrodes, 9-volt battery, hard plastic carrying case and instruction manual
- Power Source: 9-volt battery
- Channel: Dual, isolated between channels
- Pulse Width: 250 microseconds
- Pulse Rate: 5, 30, 100 Hz

### PRODUCTS

### Description

PRODUCTS	Description
DE5030	EMS 5.0 - Analog EMS 32EA/CV