

## Lead Wire Extender

HCPCS:

- The Lead Wire Extender connects to the ends of existing lead wires to add an additional 5' of length for comfortable TENS therapy at a distance.
- The extender is compatible with most standard single-pin lead wires and TENS and EMS units.
- Adds additional 5 feet to length of existing single-pin lead wires
- Compatible with most standard TENS and EMS units
- Quantity: 1 pair/pack

<b>PRODUCTS</b>	<b>Description</b>
WW3766	Lead Wire Extender 1EA/