

Relief Wrap Conductive Brace with 2" x 2" Electrodes, Large

HCPCS:

- The Relief Wrap Conductive Brace helps users receive electrode therapy in hard to reach areas without requiring the help of others.
- It can be customized to fit patient's needs based on treatment protocol, accommodating either 4 electrodes for TENS and EMS therapy or 8 electrodes for interferential therapy.
- Comes with 4-pack of attachable 2" x 2" electrodes (#EP2020RW2)
- Size large, fits approximately 30" to 40" waists; also available in small, medium, and extra-large
- Easy to maintain, hand washable
- Compatible with most devices
- To independently receive electrode therapy in even the hardest to reach areas, purchase the Relief Wrap Conductive Brace.

| PRODUCTS | Description |
|----------|------------------------------|
| GF5240 | Relief Wrap - Large 2x2 1EA/ |