

Self-Adhesive Electrodes, 1.5" x 1.5" Tan Cloth in Foil Pouch

HCPCS:

- The Tan Cloth Electrodes are made with the highest clinical standards and fit most TENS, EMS, IF, and other electro-stim applications.
- These highly conductive electrodes deliver more current to better distribute across entire pad.
- With premium tan cloth and a full-coverage wire connector, they are durable and safe for any user.
- Self-adhesive pad for ease of use
- Flexible carbon layer allows for maximum charge delivery and comfort on a variety of application sites
- Durable and comfortable tan cloth backing
- Durable and reliable full-coverage wire connector
- Sealed in a tamper-proof foil pouch
- Size: 1.5" x 1.5"
- Quantity: 4 electrodes/pack
- Can be reused up to 15 times per pack

PRODUCTS

Description

E1F1515TC2

Electrode - Foil - 1.5x1.5 - 400EA/CV
