

## Self-Adhesive Electrodes, 1.5" x 1.5" White Cloth in Poly Bag

HCPCS:

- The White Cloth Electrodes are made with the highest clinical standards and fit most TENS, EMS, IF, and other electro-stim applications.
- These highly conductive electrodes deliver more current to better distribute across entire pad.
- With premium white cloth and a full-coverage wire connector, they are durable and safe for any user.
- Self-adhesive pad for ease of use
- Flexible carbon layer allows for maximum charge delivery and comfort on a variety of application sites
- Durable and comfortable white cloth backing
- Durable and reliable full-coverage wire connector
- Made with premium gel
- Sealed in a tamper-proof poly bag
- Size: 1.5" x 1.5"
- Quantity: 4 electrodes/pack, 10 packs/poly bag

### PRODUCTS

### Description

---

E1P1515WC2	Electrode - Poly - 1.5x1.5 - 400EA/CV
------------	---------------------------------------

---