

Self-Adhesive Electrodes, 2" x 3.5" Blue Cloth with Tyco Gel, Poly Bag

HCPCS:

- The Blue Cloth Electrodes are made with the highest clinical standards and fit most TENS, EMS, IF, and other electro-stim applications.
- These highly conductive electrodes deliver more current to better distribute across entire pad.
- They are self-adhesive and easy to use, and with premium blue cloth and a full-coverage wire connector, they are durable and safe for any user.
- Self-adhesive pad for ease of use
- Flexible carbon layer allows for maximum charge delivery and comfort
- Durable and comfortable premium blue cloth backing
- Durable and reliable full-coverage wire connector
- Tyco gel is self-adhesive and made in the USA
- Size: 2" x 3.5"
- 4 electrodes/pack; 10 packs/polybag
- Can be reused up to 15 times per pack

PRODUCTS

Description

EP2035BC2	Electrode - Poly - 2x3.5 - 320EA/CV
-----------	-------------------------------------
