

Self-Adhesive Electrodes, 3" Round White Foam, Foil Pouch

HCPCS:

- The 3" Round White Foam Electrodes are made with the highest clinical standards and fit most TENS, EMS, IF, and other electro-stim applications.
- These highly conductive electrodes deliver more current to better distribute across entire pad.
- With premium white clinical foam and a full-coverage wire connector, they are durable and safe for any user.
- Self-adhesive pad for ease of use
- Flexible carbon layer allows for maximum charge delivery and comfort
- Durable and comfortable white foam backing
- Durable and reliable full-coverage wire connector
- Tyco gel is self-adhesive and made in the USA
- Sealed in a tamper-proof foil pouch
- Size: 3" Round
- Quantity: 4 electrodes/pack, 10 packages/bag

PRODUCTS

Description

EF3000WF2

ELECTRODE - FOIL - 3" ROUND 1EA/
