

## TheraSonic Ultrasound Gel, 5 Liter Container (1.3 gallon)

HCPCS:

- The TheraSonic Ultrasound Conductive Gel is acoustically correct for a broad range of ultrasound frequencies, making it compatible with any ultrasound device.
- It is recommended for all procedures in which a viscous gel is required.
- Acoustically correct for the broad range of ultrasound frequencies
- Can be used for a variety of purposes, including fetal doppler baby heart monitor, doppler ultrasound, ultrasonic cavitation gel, and therapeutic ultrasound procedures
- Spreads evenly and won't stain clothing or damage transducers
- Water-soluble, hypoallergenic, bacteriostatic, and non-irritating
- Recommended for all procedures in which a viscous gel is required
- 5-liter (1.3-gallon) container with 8 oz.
- reusable bottle for easy refills; available in a variety of other sizes
- For an effective viscous gel for ultrasound procedures, purchase the TheraSonic Ultrasound Conductive Gel.

### PRODUCTS

### Description

---

LS5255	TheraSonic Conductive Gel, 5L 1EA/
--------	------------------------------------

---