

Twin Stim TENS and EMS Combo 2nd Edition

HCPCS:

- Our most popular digital electrotherapy unit.
- The Twin Stim 2nd Edition is an exact combination of our TENS 7000 and EMS 7500, which gives the user the choice of five Modes of TENS and two Modes of EMS.
- Features include, dual channels, seven modes, timer, patient compliance meter, new stronger 100 mA output, safety amplitude cover, sturdy metal belt clip and rubber side railings for maximum handling.
- Five TENS Modes: Burst / Normal / Modulation / SD1 / SD2
- Two EMS Modes: Synchronous, Alternate
- System includes device, lead wires, four self-adhesive reusable electrodes, 9-volt battery, hard plastic carrying case and instruction manual
- Patient compliance meter can store up to 60 sets of operation manual
- Power Source: 9-volt battery
- Pulse Width: 50-300 s, adjustable
- Pulse Rate: 2-150 Hz, adjustable
- Pulse Ramp: 1-8 seconds

PRODUCTS

Description

DS2202	Twin Stim 2nd - TENS/EMS 32EA/CV
--------	----------------------------------
