

Pediatric Fingertip Pulse Oximeter

HCPCS:



FIGURE A

- Four-way OLED with six display modes
- Prescription required
- Portable, non-invasive spot check of oxygen saturation and pulse rate
- Uses two AAA batteries
- Low battery indicator
- Comes in retail box
- Easy operation auto-off feature
- Lanyard, silicon case, and hang bag included
- Child-friendly design
- Adjustable brightness levels from 1 - 10
- Finger size range: 0.3" - 1.0"

PRODUCTS	Description
18707	Pediatric Pulse Oximeter, 1/ea
SPECIFICATIONS	Description
Warranty	Limited 2 Year