

## Meridian Ultra-Care Excel System

HCPCS:



- Includes mattress, cover and pump
- Features a 5" alternating pressure upper surface combined with an integrated 3" foam base to provide a full eight inches of therapeutic support with bottom-out protection in the event of a power failure
- Separate compartment for foam base inside polyurethane cover prevents soiling
- System includes mattress, cover and pump
- Mattress and pump come with a 2-year warranty
- Quick disconnect cells for easy cell replacement
- Air cells made with nylon/polyurethane for increased durability and reliability
- CPR quick deflate located at the head of the mattress
- Quilted cover, moisture resistant, vapor permeable nylon; low shear/low friction
- Pump is an 8 liter-per-minute pump in an easy-to-use design
- 10 minute cycle time

**PRODUCTS**

**Description**

PRODUCTS	Description
MER-4500	8" APM-8LPM PUMP W/FOAM BASE 1EA/