

SATINAIR APM-5/3 W/ANALOG PUMP

HCPCS:



FIGURE A



FIGURE B



FIGURE C



FIGURE D

- The Meridian SatinAir Foam (5"air/3"foam) APM system is a complete therapy mattress system.
- Designed to provide alternating pressure and low air loss, the SatinAir Foam optimizes pressure redistribution and helps to manage skin maceration.
- This system is the perfect option where a base foam layer is critical to prevent patients bottoming out during power outages.
- The powerful 9 liter-per-minute pump offers static inflation for transfers or patient care and convenient settings by patient weight.
- The nylon zipper cover is low shear, fluid-resistant and vapor permeable.
- The emergency CPR valve allows for quick deflation if necessary.
- The Meridian SatinAir Foam offers innovative design for use in nursing homes, hospitals and the home care environment.

PRODUCTS

Description

SA-SAF350A

SATINAIR APM-5/3 W/ANALOG PUMP 1EA/