

SATINAIR APM-8"W/ANALOG PUMP

HCPCS:



FIGURE A



FIGURE B



FIGURE C



FIGURE D



FIGURE E

- Meridian Medical 8" APM System
- The Meridian SatinAir 8" APM system is a complete therapy mattress system.
- Designed to provide alternating
- pressure and low air loss, the SatinAir optimizes pressure redistribution and helps to manage skin maceration.
- The powerful 9 liter-per-minute pump offers static inflation for transfers or patient care and convenient settings by
- patient weight.
- The nylon zipper cover is low shear, fluid-resistant and vapor permeable.
- The emergency CPR valve
- allows for quick deflation if necessary.
- The Meridian SatinAir offers innovative design for use in nursing homes, hospitals and the home care environment.
- FEATURES

PRODUCTS	Description
SA-SSC350A	SATINAIR APM-8"W/ANALOG PUMP 1EA/