

## ULTRACARE BAR APM-42"W/DIGPUMP

HCPCS:



FIGURE A



FIGURE B



FIGURE C



FIGURE D

- Meridian Medical 10" Bariatric APM System
- The deluxe Meridian Medical UltraCare Bariatric 10 APM is a combination mattress system for patients up to 800 lbs.
- This durable pressure reduction mattress provides alternating pressure to help prevent skin breakdown.
- Plus, the cell-on-cell construction provides bottom-out patient protection in the event of a power failure!
- Advanced digital 12 LPM pump includes easy weight and adjustment for optimal performance and comfort.
- FEATURES
- Includes cell-on-cell mattress, cover and digital pump
- Cell-on-cell design prevents the patient from bottoming out during a power outage
- Air cells can be easily removed and replaced
- 10" air cells made with nylon/polyurethane for increased durability and reliability
- CPR quick deflate located at the head of the mattress

### PRODUCTS

### Description

UC-SCC800D

ULTRACARE BAR APM-42"W/DIGPUMP 1EA/