

ELECTRONIC APM PUMP

HCPCS:

- This Digital High-output Pump is included on select APM systems or an upgrade pump for select Ultra-Care APM's.
- It offers both alternating pressure and static operation modes, plus the high-efficiency dual compressor increases patient comfort and the systems weight capacity.
- Digital pump is high-flow, high-efficiency in an easy-to-use design
- Features dual compressor for increased patient comfort and add 100 lbs to weight capacity of select APM's
- Offers both alternating pressure or static operation modes; will automatically reset from static to alternating after 45 minutes
- Includes 10, 15, or 20-minute cycle times and five comfort settings
- Has two ports and is compatible with most Ultra-Care APM's
- For an upgrade from the standard pump, choose the Digital APM Pump.
- **WARNING:** This product can expose you to chemicals including Lead, which is known to the State of California to cause cancer or birth defects or other reproductive harm.
- For more information, visit www.P65Warnings.ca.gov .

PRODUCTS

Description

4253E

DIG. APM PUMP, DUAL COMPRESSOR 4EA/CV