

MAGICAIR PUMP (FOR APM'S)

HCPCS:

- The quick-fill MagicAir Pump is designed to quickly inflate and deflate Meridian Low Air Loss Mattresses reducing your patient set-up time by 65%, and allowing you time to attend to your patient's other needs.
- This fast 51 CFM flow pump allows your technician to get the mattress up and ready for therapeutic use in minutes.
- Lightweight, compact and with steel construction - assembled in the USA
- Pumps air at a rate of 51 CFM
- Quieter as compared to standard plastic air pumps
- Weight: 3 lbs.
- Dimensions: 13" length x 5.5" wide height
- 12 foot cord
- One year warranty
- Replaces item 90501
- Don't forget the Quick Fill Valve 90153 (sold separately) - Get Meridian APM's patient-ready in minutes!

PRODUCTS

Description

109-516811	MAGICAIR® PUMP (FOR APM'S) 6EA/CV
------------	-----------------------------------