

Replacement Ventilated 8" Cell for APM

HCPCS:

- This is the replacement ventilated 8" cell for Roscoe 8" Alternating Pressure Mattresses.
- It is easy to install and connect to mattress to replace cells that may be worn out or damaged.
- Replaces APM cells that may be worn out or damaged
- Constructed of durable nylon/polyurethane
- Easy to remove and replace; easily snaps into and connects to mattress
- Size: 8"
- Ventilated; non-ventilated cell also available
- For use with Roscoe Alternating Pressure Mattresses (Genesis III, Altus, etc.)
- Sold individually

PRODUCTS	Description
90145	Replacement 8 Vented Cell 1EA/