

Foam Ring Cushion

HCPCS:



FIGURE A

- Provides excellent comfort and support while seated
- Open center is designed to provide maximum pressure relief
- Comes with durable, washable, zippered jacquard cover
- Made from soft, premium quality therapeutic foam
- Retail Packaged
- Dimensions: 16" (L) x 12 75" (W) x 2 75" (H)
- Weight Capacity: 250lbs

PRODUCTS	Description
RTLPC23388	Invalid Ring - 16",1EA,6/CS
SPECIFICATIONS	Description
Material	Jacquard Knit Cover