

Meridian Optimum Pommel Cushion (18" x 16" x 3")

HCPCS:

- The Meridian Optimum Pommel Cushion provides additional comfort and pressure relief to patients spending long periods of time in their wheelchair.
- It features a contoured base and pommel, with four grades of foam for optimum support.
- Features contoured base and pommel
- Four grades of foam for optimum support and comfort
- Inner core has a liner to protect foam
- Cover is a 4-way stretch polyester/polyurethane
- Moisture repellent and low shear/friction
- Outer cover color may be blue/black or all black
- Dimensions: 18" x 16" x 3"
- Weight capacity: 275 lbs
- One-year warranty

PRODUCTS

Description

HDC-POM1816	POMMEL CUSHION, 18X16 4EA/CV
-------------	------------------------------