

Molded Foam General Use Wheelchair Cushion

HCPCS: E2601-1,E2601-1



FIGURE A

- A unique construction of high-resilient molded foam provides optimal support, stabilization and pressure redistribution
- Designed for low risk patients seeking comfort and support
- Removable, zippered, low-shear stretch cover is vapor permeable and water-resistant
- Cover has a "Non-Slip" vinyl bottom for comfort and safety
- Convenient 12" tie straps secure the cushion to wheelchair
- Meets flammability standard TB117

PRODUCTS	Description
14880	W/C Cushion Molded,16x16x2, 1/ea, 6/cv
14881	W/C Cushion Molded, 20x16x2, 1/ea, 6/cv
14887	W/C Cushion Molded, 18x16x2, 1/cv, 6/ct
14908	W/C Cushion Molded, 18x18x2, 1/ea, 6/cv
14909	W/C Cushion Molded, 20x18, 1/ea, 6/cv
SPECIFICATIONS	Description
Warranty	18 Months Limited
Weight Capacity	250 lbs 300 lbs
Cover Material	Fluid-Resistant Stretch Nylon
CA Prop 65	[WARNING: This product can expose you to chemicals that are known to the State of California to cause cancer and birth defects or other reproductive harm. For more information go to www.P65Warnings.ca.gov .]