

POMMEL CUSHION, 20X16

HCPCS:

- The Meridian Optimum Pommel Cushion provides additional comfort and pressure relief to patients spending long periods of time in their wheelchair.
- It features a contoured base and pommel, with four grades of foam for optimum support.
- Features contoured base and pommel
- Four grades of foam for optimum support and comfort
- Inner core has a liner to protect foam
- Cover is a 4-way stretch polyester/polyurethane
- Moisture repellent and low shear/friction
- Dimensions: 18" x 16" x 3"
- Weight capacity: 275 lbs
- One-year warranty
- Purchase your Meridian Optimum Pommel Cushion for added comfort and pressure relief on your wheelchair seat.

| PRODUCTS | Description |
|-----------------|------------------------------|
| HDC-POM2016 | POMMEL CUSHION, 20X16 4EA/CV |