

Reusable Nebulizer Kit, 10 kits / case

HCPCS: A7005-1,A7005-1



FIGURE A

- Anti-spill "T" design
- Leak-proof and secure seal with simple 1/4 turn assembly design
- Very low residual volume (0.7cc)
- High medication bowl capacity (10cc)
- Graduated volume (cc) markings on bowl
- Nebulizes in upright positions or at angles up to 45°
- Dishwasher safe

PRODUCTS

Description

| PRODUCTS | Description |
|----------|--------------------------|
| 18254 | Reusable Neb Kit 10/case |