

## Airial Holding Chamber

HCPCS: A4627-1,A4627-1



FIGURE A

- Momentarily holds aerosol medication in chamber until patient is ready to inhale
- Accommodates most CFC and HFA meter dose inhalers
- Helps to eliminate waste of medication
- Minimizes need for hand-breath coordination
- Maximizes drug delivery
- Easy-to-use and clean
- Anti-static properties and Latex Free

PRODUCTS	Description
MQ8000	Airial chamber 1/ea
MQ8000L	Airial Chamber for MDI w/ Lg Mask 1/ea
MQ8000M	Airial Chamber for MDI w Med Mask 1/ea
MQ8000S	Airial Chamber for MDI w Sm Mask 1/ea