

PulmoMate® Compressor Nebulizer System

HCPCS: E0570



- The sleek case design of the PulmoMate® Compressor Nebulizer System fits easily into luggage, backpacks and diaper bags for easy transport The DeVilbiss PulmoMate features an updated compressor/motor for long-term durability and performance
- Robust quality and durability
- Consistent performance over the life of the unit
- Metal frame motor construction
- Integrated handle for easy carrying
- Large built-in compartment for storing nebulizers and tubing
- Intuitive design is easy-to-use and maintain

PRODUCTS	Description
4650D	PulmoMate w/ Disp. Nebulizer
SPECIFICATIONS	Description
Warranty	Limited 5 year
Dimensions	7.0" (W) x 3.8" (H) x 13.0" (D)
Free Air Flow	13 lpm Typical
Weight	5.3 lbs
ACCESSORIES	Description
3655D-621	VixOne Disp. Nebulizer, CS 50, CS
DL1050	PEDIATRIC CHARACTER MASK (50-PK), CS
5650D-602	FILTER PKG - MM, 5/PK, PKG
4650D-605	STORAGE COMPARTMENT DOOR ASSY, EA
5650D-612	FILTER KIT(CAP&FILTER) (6/PK), PKG
4650D-609	PUSH-IN FEET (4-PK), PKG
5650D-601	FILTER PKG - DFT (25-PK), PKG
4650D-611	TUBE CONNECTOR W/PAL NUT 6-PKG, EA