

Pacifica Elite Nebulizer

HCPCS: E0570-1,E0570-1



FIGURE A



FIGURE B



FIGURE C

- Retail-packaged compressor neb is a workhorse with higher liter flows, shorter treatment times and quieter operation than similar economy compressor nebs
- Ultra compact design makes it extremely portable
- Powerful piston pump is efficient, dependable and quiet
- Delivers consistent particle size of <5 microns and fast treatment times
- Retail packaged
- CAUTION: Federal law restricts this device to sale by, or on the order of a licensed healthcare practitioner.

PRODUCTS	Description
18070	Pacifica Elite Nebulizer,6RTL/CS, 6/cv
18071	Pacificaelite w/Reuse&DisposeNebKit,6/CS
SPECIFICATIONS	Description
Warranty	Limited 5 Year
Electrical Requirements	115 VAC 60Hz
Liter Flow	8 lpm
Method of Nebulization	Piston Compressor
Operating Pressure	10 psi
Partical Size	0.5 micron to 5 micron
Weight	4.1 lbs.
ACCESSORIES	Description
18031	Nebulizer Carry Bag, 1/ea, EA
18254	Reusable Neb Kit 10/case, CS
MASK 001A	Adult Aerosol mask CS/50, CS
NEB KIT 500	Neb Kit,Inc Mouth pc&Tub50/cs, CS
18090F	NebFilter 10/bag, BAG