

Power Neb Ultra Nebulizer

HCPCS: E0570, A7003, A7005, A7013



- A competitive edge price point for providers as one of the smallest compressor nebs at a mere 3.2 lbs
- Built-in neb kit holder for hands-free use
- Delivers consistent particle size of <5 microns and fast treatment times
- Complete with disposable neb kit, and optional reusable neb kit or pediatric aerosol mask
- Powerful piston pump is efficient, dependable and quiet
- CAUTION: Federal law restricts this device to sale by, or on the order of a licensed healthcare practitioner.

PRODUCTS	Description
18080	Power Nebulizer, Ultra, 1/ea, 6/cv
18081	PowerNebUltra/Reuse&DispNeb, 1/ea, 6/cv
18082	PowerNeb, Ultra w/2 Reusable kits,6/cs
18260-PDMASK	Ped Mask for 18260 and Adapter, 1/ea

SPECIFICATIONS	Description
Warranty	Limited 5 Year
Liter Flow	, 8 lpm
Method of Nebulization	Piston Compressor
Operating Pressure	, 10 psi
Particle Size	0.5 µm to 5 µm
Weight	3.2 lbs.

ACCESSORIES	Description
18031	Nebulizer Carry Bag, 1/ea, EA
18254	Reusable Neb Kit 10/case, CS
NEB KIT 600	NebKit,Inc Adlt Mask,Tub 50/cs, CS
MASK 001A	Adult Aerosol mask CS/50, CS
MASK 001P	Pediatric Aerosol mask CS/50, CS
NEB KIT 700	NebKit,IncPediMask,Tub ,50/cs, CS
NEB KIT 500	Neb Kit,Inc Mouth pc&Tub50/cs, CS
18260	JetNeb Plus Breath Enhanced Neb,12/CS, CS
18080F	18080.18081.18082, NebFilter, 10/bag, BAG