

Model 15-RD Atomizer

HCPCS:



- First produced in the early 1900s, features the same adjustable swivel tip and vented nasal guard that made it as popular then as it is today
- Metal atomizer with adjustable tip for directing spray up or down to reach inaccessible areas of the throat
- Can be used with or without nasal guard Used for nose and throat treatment hand bulb
- Basically, all types of oil or water solutions can be used that are compatible with nickel plating
- Used in colleges and laboratories for experimental testing
- Metal banded tubes and spray tips (sold separately) are used by research laboratories

PRODUCTS	Description
15-RD	Glass Atomizer w/ Plastisol Bulb

ACCESSORIES	Description
115-BT-453-605	BRASS TUBING ASSEMBLY 5 PACK, EA
114-S-605	WASHER HERCULES 10-481 5-PK, PKG
115-QB	1-OZ CLEAR GLASS, EA
115-CLN-605	POINT 5-PK, PKG
127-G-605	2" POLY OUTLET TUBE 5-PK, PKG
12-RLN-605	NASAL GUARD 5PK, PKG
1-RD-605	PLASTISOL BULB (5/32" OPENING) 5 PK SS, PKG
1-RD	PLASTISOL BULB (5/32" OPENING), EA