

## Model 163 Atomizer

HCPCS:



- The Model 163 Atomizer offers one of our most consistent sprays for aqueous or viscous solutions. The adjustable tip helps direct the fine spray precisely where needed, and the all-metal top is autoclavable.
- Designed for quick change and replacement of spray tube.
- Tight-fitting "O" ring ensures consistent spray.
- Metal spray tube may be sterilized after each use.
- All types of oil or water solutions can be used that are compatible with rhodium and nickel plating.
- Adjustable tip permits up or down spray for aqueous solutions of heavier viscosity.
- May be operated by a No 630 Bleeder-Type cut-off assembly or a No 1 bulb.
- Used for decubitus ulcer treatments, debridement of wounds, and wound cleansing.
- Used by E.N.T doctors, research laboratories, veterinarians, burn centers, wound care facilities, emergency rooms.

### PRODUCTS

#### Description

PRODUCTS	Description
163-RD	Glass Atomizer w/ Metal Top

### ACCESSORIES

#### Description

ACCESSORIES	Description
150-QB	1.4 OZ CLEAR GLASS, EA
163-S-605	O- RING PARKER SIZE 2-7, 5 PK, PKG
163-T-605	O-RING SIZE 2-9, 5 PK, PKG
115-CLN-605	POINT 5-PK, PKG
251-S-605	WSHER 10-481 STK EQ 5 PK SEPE SALE PKG, PKG
150-OT-605	OUTLET TUBE ASSY 5 PK SEPE SALE PKG, PKG
630-605	CUTOFF WITH TUBING 5-PK, PKG